

IHE Integration Statement			
Vendor	Product Name	Version	Date
CGM Clinical Österreich	Integration Package (i.Pack)	4.1.2.x	3/17/21

This product implements all transactions required in the IHE Technical Framework to support the IHE Integration Profiles, Actors and Options listed below :

Integration Profile	Actor	Integration Profile Option
Audit Trail and Node Authentication	Secure Application	None
Audit Trail and Node Authentication	Secure Application	ATX: TLS Syslog
Audit Trail and Node Authentication	Secure Application	STX:TLS 1.2 Floor using BCP195
Consistent Time	Time Client	None
Community Dispense	Content Consumer	View Option
Community Dispense	Content Consumer	None
Document Metadata Subscription	Document Metadata Notification Recipient	None
Document Metadata Subscription	Document Metadata Subscriber	None
Document Metadata Subscription	Document Metadata Subscriber	Patient-Independent Subscription
Community Pharmaceutical Advice	Content Consumer	View Option
Community Pharmaceutical Advice	Content Consumer	None
Patient Demographics Query	Patient Demographics Consumer	None
Patient Demographics Query for Mobile	Patient Demographics Consumer	None
Patient Demographic Query HL7 V3	Patient Demographics Consumer	None

Internet address for vendor's IHE information :

[http://www.cgm.com/at/krankenhaus\\_klinik\\_1/kommunikationsloesungen\\_1/cgm\\_ihe\\_1/downloads\\_8/downloads\\_ihe.de.jsp](http://www.cgm.com/at/krankenhaus_klinik_1/kommunikationsloesungen_1/cgm_ihe_1/downloads_8/downloads_ihe.de.jsp)

Links to Standards Conformance Statements for the Implementation

<b>HL7</b>	
<b>DICOM</b>	

Link to general information on IHE : [www.ihe.net](http://www.ihe.net)

IHE Integration Statement			
Vendor	Product Name	Version	Date
CGM Clinical Österreich	Integration Package (i.Pack)	4.1.2.x	3/17/21

This product implements all transactions required in the IHE Technical Framework to support the IHE Integration Profiles, Actors and Options listed below :

Integration Profile	Actor	Integration Profile Option
Patient Identifier Cross-referencing for MPI	Patient Identity Source	None
Patient Identifier Cross-referencing for MPI	Patient Identifier Cross-reference Consumer	None
Patient Identifier Cross-reference for Mobile	Patient Identifier Cross-reference Consumer	None
Patient Identifier Cross-Reference HL7 V3	Patient Identity Source	None
Patient Identifier Cross-Reference HL7 V3	Patient Identifier Cross-reference Consumer	None
Community Prescription	Content Consumer	View Option
Community Prescription	Content Consumer	None
Cross-Enterprise Document Reliable Interchange	Document Source	None
Cross-Enterprise Document Reliable Interchange	Metadata-Limited Document Source	None
Cross-Enterprise Document Sharing	Document Administrator	None
Cross-Enterprise Document Sharing	Document Consumer	None

Internet address for vendor's IHE information :

[http://www.cgm.com/at/krankenhaus\\_klinik\\_1/kommunikationsloesungen\\_1/cgm\\_ihe\\_1/downloads\\_8/downloads\\_ihe.de.jsp](http://www.cgm.com/at/krankenhaus_klinik_1/kommunikationsloesungen_1/cgm_ihe_1/downloads_8/downloads_ihe.de.jsp)

Links to Standards Conformance Statements for the Implementation

<b>HL7</b>	
<b>DICOM</b>	

Link to general information on IHE : [www.ihe.net](http://www.ihe.net)

IHE Integration Statement			
Vendor	Product Name	Version	Date
CGM Clinical Österreich	Integration Package (i.Pack)	4.1.2.x	3/17/21

This product implements all transactions required in the IHE Technical Framework to support the IHE Integration Profiles, Actors and Options listed below :

Integration Profile	Actor	Integration Profile Option
Cross-Enterprise Document Sharing	Document Source	None
Cross-Enterprise Document Sharing	Document Source	Folder Management
Cross-Enterprise Document Sharing	Document Source	Document Addendum
Cross-Enterprise Document Sharing	Document Source	Document Transformation
Cross-Enterprise Document Sharing	Document Source	Document Replacement
Cross-Enterprise Document Sharing for Imaging	Document Consumer	None
Cross-Enterprise Document Sharing for Imaging	Imaging Document Consumer	None
Cross-Enterprise Document Sharing of Scanned Documents	Content Consumer	View Option
Cross-Enterprise Document Sharing of Scanned Documents	Content Consumer	None
Cross Enterprise User Assertion	X-Service User	None
Cross Enterprise User Assertion	X-Service User	Subject-Role
Cross Enterprise User Assertion	X-Service User	PurposeOfUse

Internet address for vendor's IHE information :

[http://www.cgm.com/at/krankenhaus\\_klinik\\_1/kommunikationsloesungen\\_1/cgm\\_ihe\\_1/downloads\\_8/downloads\\_ihe.de.jsp](http://www.cgm.com/at/krankenhaus_klinik_1/kommunikationsloesungen_1/cgm_ihe_1/downloads_8/downloads_ihe.de.jsp)

Links to Standards Conformance Statements for the Implementation

<b>HL7</b>	
<b>DICOM</b>	

Link to general information on IHE : [www.ihe.net](http://www.ihe.net)