


Control .C1 Low Chemistry

Test: Glucose (mg/dL)



Lot #: CHEM005

Date/Time	Result	Tech	Flag	Comments
Expected Mean: 65.0		2SD: 6.0		Expires: 6/19/2013
9/1/2012 10:14 AM	67.0	BR		
9/2/2012 10:00 AM	68.0	BR		
9/3/2012 10:09 AM	66.0	BR		
9/4/2012 10:18 AM	62.0	BR		
9/5/2012 10:27 AM	66.0	BR		
9/6/2012 9:49 AM	67.0	BR		
9/7/2012 9:58 AM	64.0	BR		
9/8/2012 10:06 AM	68.0	BR		
9/9/2012 10:17 AM	68.0	BR		
9/10/2012 10:25 AM	64.0	BR		
9/11/2012 10:35 AM	62.0	BR		
9/12/2012 10:03 AM	67.0	BR		
9/13/2012 10:15 AM	62.0	BR		
9/14/2012 10:08 AM	61.0	BR		
9/15/2012 10:18 AM	55.0	BR	> 3SD	recalibrated- lab
9/16/2012 10:26 AM	67.0	BR		
9/17/2012 10:01 AM	62.0	BR		
9/18/2012 10:09 AM	68.0	BR		
9/20/2012 9:37 AM	66.0	BR		
9/21/2012 10:00 AM	66.0	BR		
9/22/2012 10:07 AM	62.0	BR		
9/23/2012 9:24 AM	66.0	BR		
9/24/2012 9:44 AM	67.0	BR		
9/25/2012 9:18 AM	63.0	BR		

Low Chemistry **Calculated Mean:** 64.8 **2SD:** 6.3 **%CV:** 5 **n:** 24 **Reviewed By:** _____